

KC DATA – Automatization collection and analysis

User Manual

Updated 2021-01-19


Contents

1. Introduction	3
2. Interface guide	4
Start Window	4
Action 1: Check data	4
Action 2: Add a new Project.....	4
Action 3: Open Project	6
Action 4: Delete Project	6
Project Window	8
Action 1: New Season	9
Action 2: View info of a Season.....	10
Action 3: New experiment	10
Action 4: View Data from an experiment	11
Action 5: Delete an Experiment	12
Action 6: Analyze data	12
Experiment Window	13
Action 1: Modify Experiment Attributes.....	13
Action 2: Open a Point	14
Action 3: Add a Point	14
Action 4: View Data.....	14
Action 5: Delete Data	14
Action 6: Analyze data	14
Action 7: Add data.....	15
DATA TYPES.....	19
SCADA	19
GC data.....	20
SPA data	20
Other Files	21
Point Window	22
Action 1: Modify info	22
Action 2: View Data.....	22
Action 3: Delete Data	22
Action 4: Analyze data	23
Action 5: Link data.....	23
3. How to implement code changes	24

1. Introduction

This program is used as a server to collect all the data gathered in the *Kraftcentralen (KC or Power Central)* during research. Its purpose is to store data and facilitate its analysis.

The interface is made of 4 main windows: Start, Project, Experiment and Point. The structure is as shown in the figure below. A project, for example investigating the thermochemical process for Polyethylene, consist of various seasons (ex. winter of 2020-2021), which can have one or more experiments (commonly 1 day). One experiments could be for example investigating the effect of temperature and can contain one or more points of approx. 1-hour duration, for instance two different temperatures.




As January 19th, 2021, all data collection is finished and automatized for SCADA type (i.e trend obtained from the operation conditions in the power plant), GC type (for the gas analyzed in Gas chromatography, both for the G101 and for inferno in the gasifier, for dry gas and for the mass balance, respectively) and for SPA type (i.e. data coming from the GC-430 including mainly aromatics compounds), while the analysis will be implemented later.

2. Interface guide

Start Window

The start window includes a summary table of all the data available divided by its project.



The screenshot shows a window titled "Projects KC" with standard window controls (minimize, maximize, close). The main content area has a blue header with the text "KC - data". Below the header is a table with the following data:

	Index	Project Name	Timeframe	Description	Resp
1	0	Proj0		this is project 0	resp0
2	1	Proj1		this is project 1	resp1
3	2	Proj2		this is project 2	resp2

Below the table is a horizontal scrollbar. At the bottom of the window are four buttons: "DELETE PROJECT", "CHECK DATA", "ADD PROJECT", and "OPEN PROJECT".


Action 1: Check data

If you want to check, open or download raw data click *Check data*

This function is under construction 

Action 2: Add a new Project

If you need to add a project click *Add Project*:



New project


Project Name

Description

Responsible

Create Cancel

Type the project name, description and responsible and click *Create*:



New project

Project Name

Description

Responsible

Create Cancel


The new project line will now appear in the Start Window:

Project 1		user guide ...	me
-----------	--	----------------	----

If you don't want to create a new project click *Cancel*.

Action 3: Open Project

If you want to see or modify the information of a project, select the project and click *Open Project*:




The Project Window will open (see Project Window).

Action 4: Delete Project


If you want to delete a Project, click the experiment and click *Delete Exp.*

Warning! It is not recommended to delete an experiment; the analysis may also be deleted.

An are you sure message will appear; the project will be deleted if you press yes.




If no project is selected an error message will appear.




Project Window

The window includes the main information of the project: name and description; and a list with all the seasons and experiments available; and allows to do multiple actions.




To deploy more information about the experiments and points click the arrows as shown in the figure below:



Action 1: New Season


If you need to add a new Season click *New Season* and add name and decryption and click *Create*:



The new Season line will now appear in the Project Window.


If you don't want to create a new project click *Cancel*.

Errors will appear if the information needed is not available.



Action 2: View info of a Season

If you want to see the existing information of a Season Click *View Info*:



View info Season X

From: 1/15/2021 To: 1/20/2021

Description: This is a test Season

Experiments list:

- Exp A: high temperature
- Exp B: Low temperature
- Exp C: Circulation
- Exp D: asadfad


Buttons: Cancel, OK

You can change the description only and click *ok*, or just see the information available and close the window or click *Cancel*.

The date is automatically filled with the minimum and maximum dates of the experiments.

Action 3: New experiment

If you want to add a new experiment, click the season where the experiment takes place and click *New Experiment*:



NAME	Date Ini	Date End	Ent.	Fuel Type
✓ Ses_p1_s1	2020-10-02 10:00:00	2020-10-18 12:00:00	1	Polyethylen
> exp0	2020-10-2 10:00:00	2020-10-18 12:00:00	4	Polyethylene
✓ Ses_p1_s2	2020-10-08 10:00:00	2020-10-18 12:00:00	1	PVC
> exp0	2020-10-8 10:00:00	2020-10-18 12:00:00	3	PVC
✓ Ses_p1_s3	2020-10-01 10:00:00	2020-10-28 12:00:00	4	Mix
> exp0	2020-10-1 10:00:00	2020-10-18 12:00:00	1	PVC
> exp1	2020-10-10 10:00:00	2020-10-28 12:00:00	3	Textiles
> exp2	2020-10-6 10:00:00	2020-10-21 12:00:00	3	Polyethylene
> exp3	2020-10-2 10:00:00	2020-10-20 12:00:00	2	Polyethylene
Season 2...	ND	ND	0	ND

Buttons: NEW SEASON, VIEW INFO, NEW EXPERIMENT, VIEW DATA, DELETE EXP, ANALYZE DATA

A new window will open with the information needed for the experiment, as shown below.

The season selected will appear, but another season within the project can be selected too. If you need a new season you can click *New Season* to create a new one (See action 1 New season).

The experiment numbering is automatic, and it also take the information of the last experiments (bed material and fuel). Change those parameters if needed and add the date and times of the experiment. When finish Click *Create*.

It is important to comment what the experiment is about to make easier to everyone the data analysis.

If you don't want to add a new experiment, click *Cancel*.

Action 4: View Data from an experiment

If you want to open an experiment, click the experiment and click *View data*:

NAME	Date Ini	Date End	Ent.	Fuel Type
✓ Ses_p1_s2	2020-10-08 10:00:00	2020-10-18 12:00:00	1	PVC
> exp0	2020-10-8 10:00:00	2020-10-18 12:00:00	3	PVC
✓ Ses_p1_s3	2020-10-01 10:00:00	2020-10-28 12:00:00	4	Mix
> exp0	2020-10-1 10:00:00	2020-10-18 12:00:00	1	PVC
Point0	ND	ND	-	PVC
> exp1	2020-10-10 10:00:00	2020-10-28 12:00:00	3	Textiles
Point0	ND	ND	-	Textiles
Point1	ND	ND	-	Textiles
Point2	ND	ND	-	Textiles
> exp2	2020-10-6 10:00:00	2020-10-21 12:00:00	3	Polyethylene
> exp3	2020-10-2 10:00:00	2020-10-20 12:00:00	2	Polyethylene

For more information see [Experiment window](#).


Action 5: Delete an Experiment

If you want to delete an experiment, click the experiment and click *Delete Exp.*

Warning! It is not recommended to delete an experiment; the data may also be deleted.

Action 6: Analyze data

If you want to check, open or download raw data from more than one experiments click: *Analyze data.*

This function is under construction 

Experiment Window

The *Experiment window* includes the main information of the Experiment: name, start and end date and time and fuel, bed material and comments. This information can be changed if needed in *Modify Experiments attributes*. Also includes a status information (highlighted).

It also includes two tables: Points available and Data available. Including multiple actions.

The screenshot shows the 'Experiment Window' interface. At the top, it displays 'SEASON: Ses_p1_s0/ EXPERIMENT: exp0'. The main section contains fields for experiment details: Exp. Name (exp0), Date start (2020-10-2), Date end (2020-10-28), Fuel (PVC), Bed Type (silica sand), and Comments (the bed was with iron, and it was necessary to verify potential leakages). A 'Modify Exp. Attributes' button is next to these fields. To the right, there's a 'Data Loaded: 0' section with a table for 'Data available' and buttons for 'ADD DATA', 'VIEW DATA', 'DELETE DATA', and 'ANALYZE DATA'. Below this is a 'Status' section. At the bottom, there's a 'POINTS AVAILABLE: 5' section with a table for 'Points already created and/or analyzed' and buttons for 'OPEN POINT' and 'DELETE POINT'. The interface is annotated with three callouts: 'Experiment information' pointing to the experiment details fields, 'Points already created and/or analyzed' pointing to the points table, and 'Data available in the given project and given season' pointing to the data table.

Experiment information

Points already created and/or analyzed

Data available in the given project and given season

Action 1: Modify Experiment Attributes

In order to modify the information of an Experiment (names, dates, fuel, bed type and comments) click *Modify Exp Attributes*. The text will turn white (see image) and you can modify any attribute, when finish, click *Modify Exp Attributes* again.

This screenshot shows the 'Experiment Window' interface with the 'Modify Exp. Attributes' button highlighted. The button is located next to the experiment details fields. The text in the fields is white, indicating that the 'Modify Exp. Attributes' mode is active. The interface is annotated with three callouts: 'Experiment information' pointing to the experiment details fields, 'Points already created and/or analyzed' pointing to the points table, and 'Data available in the given project and given season' pointing to the data table.

Experiment information

Points already created and/or analyzed

Data available in the given project and given season

Action 2: Open a Point

If you want to open Point, click the Point inside the points available table and click *Open Point*.

See more information at Point Window.

Action 3: Add a Point

If you want to add a new Point, click *Add Point*:

Season 1/ Experiment 1/Point 1

Name: Point0

Comments: this is the point 0

Status: Modifying Attributes...

Data Linked

Index	Type	StartDate	EndDate	Delay	Ientry
-------	------	-----------	---------	-------	--------

VIEW DATA ANALYZE DATA DELETE DATA

Date (YYYY-MM-DD) Time (HH:MM:SS)

Date start: 2020-10-1 8:00

Date end: 2020-10-1 17:00

Link data to point

Select type

- AUTOMATIC
- SCADA
- GC1
- INFERNO
- GDA

LINK DATA

OK

This action will open an empty Point Window. Fill the information and click *Create Point* when finish. See more information at Point Window.

Action 4: View Data

The table Data available shows the existing data in the given project and given season. If you want to see the raw data click *View Data*.

Action 5: Delete Data

The table Data available shows the existing data in the given project and given season. If you want to delete the raw data click *Delete Data*.

Warning! It is not recommended to delete data the data may also be deleted. It is preferable to upload the data again, a different name will be given automatically.


Action 6: Analyze data

The table Data available shows the existing data in the given project and given season. If you want to do calculation or analysis on this data click *Analyze Data*.

This function is under construction 


Action 7: Add data

If you want to add new data (which does not exist in the table Data available), select the type of data in upper right corner and click *Add data*.




Adding SCADA & GC data


When adding SCADA search your file in your computer and click *Add Data*. Add any comment if needed (recommended to write down if any anomaly happened during the experiments). See steps below:



Add SCADA window




Search for a file



Example searching for SCADA file (i.e. YYYYMMDD trend.xls)

If the file data does not correspond with the dates included in the season an error will occur:



If the dates are correct, the data will be added successfully:

OPEN EXPERIMENT >>>> SEASON: Ses_p0_s0/ EXPERIMENT: Exp 1

1

SEASON: Ses_p0_s0/ EXPERIMENT:

Exp. Name: Exp 1 Modify Exp. Attributes

Date start: 2019-02-01 08:00:00

Date end: 2019-02-01 17:00:00

Fuel: Polyethylene

Bed Type: Olevine

Comments: this was the first experiment

Data Loaded: 0

ADD DATA

SCADA
GC1
INFERNO

	DB Type	Date Start	Date End	Delay	Co
1	SCADA_0	2019-02-01 ...	2019-02-01 ...	00:00:00	This is SCAD

< >

VIEW DATA DELETE DATA ANALYZE DATA

For GC data is the same procedure, but the delay of the GC computer (not connected to internet) with respect the SCADA computer must be added (For example if GC time is 10:05 and SCADA 10:02, the delay is +3).

Add Data (GC gas /GC INFERNO)

Add Data (GC gas /GC INFERNO)

HH:MM:SS

Time Delay (+/- from SCADA): 00:03:00

File: datafile data file path Search

Comments: This is GC_0

CANCEL ADD DATA

Adding SPA data

When adding a file type of SPA data, it will automatically fill a table with the different samples and repetitions.

List of Sheets in SPA excel file

The date corresponds with the SPA collection not the characterization

Data collected: 2019-02-01 [Change Date](#)

GPX samples filled: 0

	Sheets	(YYYY-MM-DD)	(HH:MM:SS)	Bial/Repet.	G Label
7	430_190206_190201_G21_C_46_3	2019-02-01		C_3	G21
8	430_190206_190201_G21_C_46_2	2019-02-01		C_2	G21
9	430_190206_190201_G21_C_46_1	2019-02-01		C_1	G21
10	430_190206_190201_G21_C2_45_3	2019-02-01		C2_3	G21
11	430_190206_190201_G21_C2_45_2	2019-02-01		C2_2	G21
12	430_190206_190201_G21_C2_45_1	2019-02-01		C2_1	G21

[Read Table](#) [Cancel](#) [Fill Fields](#)

Add the time of sample collection where it corresponds. You can change the dates of all the points simultaneously in the upper right corner (*Change Date*)

List of Sheets in SPA excel file

The date corresponds with the SPA collection not the characterization

Data collected: 2019-02-01 [Change Date](#)

GPX samples filled: 3

	Sheets	(YYYY-MM-DD)	(HH:MM:SS)	Bial/Repet.	G Label
26	430_190206_190201_G13_C_40_2	2019-02-01	12:10:00	C_2	G13
27	430_190206_190201_G13_C_40_1	2019-02-01	12:10:00	C_1	G13
28	430_190206_190201_G13_C2_39_3	2019-02-01	12:10:00	C2_3	G13
29	430_190206_190201_G13_C2_39_2	2019-02-01	12:10:00	C2_2	G13
30	430_190206_190201_G13_C2_39_1	2019-02-01	12:10:00	C2_1	G13
31	430_190206_190201_G12_F_38_3	2019-02-01	12:05:00	F_3	G12

[Read Table](#) [Cancel](#) [Fill Fields](#)

When clicking *Read Table*, the data will be added to the experiments:

Data Loaded: 1

[ADD DATA](#)

Select Data:
GC1
INFERNO
SPA

Index	DB Type	Date Start	Date End	Delay	This is S
0	SPA_0	2019-02-01 12:00:00	2019-02-01 12:10:00	00:03:00	This is S


[VIEW DATA](#) [DELETE DATA](#) [ANALYZE DATA](#)

Status: Ready!

DATA TYPES

Important!!! The data extracted from the KC measurement equipment should stay intact. Any modifications of variable or sheets names will interfere with the automatic analyses of the program.

SCADA



Definition: Power central process conditions

Format: Excel, two sheets: channel and data

Example channel sheet:

A	B	C	D	E	F	G	H	I
nBlock	sPrefix	nChannelNumber	sSource	sTag	sDescription	nType	sUnits	fMin
	4 C	101	CALCULATED CHANNELS	LV 210	TOTAL AIR FLOW		2 kg/s	
	4 C	102	CALCULATED CHANNELS	LV 204	PRIMARY AIR FLOW		2 kg/s	
	4 C	103	CALCULATED CHANNELS	RV 205	RECIRC. FLUE GAS FLOW		2 kg/s	
	4 C	104	CALCULATED CHANNELS	LI 203	TOT. SEC. AIR FLOW		2 kg/s	
	4 C	110	CALCULATED CHANNELS	LV 286	AIR/FLUE GAS TO PART. SEAL FLOW		2 kg/s	
	4 C	115	CALCULATED CHANNELS	LI 226	FUEL SPREADER FLOW		2 kg/s	
	4 C	120	CALCULATED CHANNELS	DD89+90	MEAN TEMP. IN BOTTOM BED		2 °C	
	4 C	121	CALCULATED CHANNELS	DC27+DC28	MEAN TEMP. IN TOP		2 °C	
0	4 C	134	CALCULATED CHANNELS	C102+C103	Gas velocity in primary zone		2 m/s	
1	4 C	135	CALCULATED CHANNELS	C101	luftshastighet i toppen		2 m/s	
2	4 C	136	CALCULATED CHANNELS	C103	Rökgasåterföringshastighet i top		2 m/s	
3	4 C	137	CALCULATED CHANNELS	C101+C103	Gas velocity in top		2 m/s	
4	4 C	153	CALCULATED CHANNELS	BUNKER NR 4	Fuel flow to gasifier		2 kg/h	
5	4 C	154	CALCULATED CHANNELS	AI-302	Mass flow of steam to gasifier		2 kg/h	
6	4 C	155	CALCULATED CHANNELS	AI-305	Mass flow of steam to seal 1		2 kg/h	
7	4 C	156	CALCULATED CHANNELS	AI 306	Mass flow of steam to seal 2		2 kg/h	
8	4 C	157	CALCULATED CHANNELS	AI 302	Density of steam (gas law)		2 kgm3	
9	4 C	158	CALCULATED CHANNELS	AI-320	Mass flow of steam to dividers.		2 kg/h	
0	1 DC	220	DATASCAN_1	26-27	DP 4380 _ 5385 mm		0 kPa	
1	1 DC	221	DATASCAN_1	27-28	DP 5385 _ 6540 mm		0 kPa	
2	1 DC	222	DATASCAN_1	28-29	DP 6540 _ 7040 mm		0 kPa	
3	1 DC	223	DATASCAN_1	29-30	DP 7040 _ 7540 mm		0 kPa	
4	1 DC	224	DATASCAN_1	30-31	DP 7540 _ 7940 mm		0 kPa	
5	1 DC	225	DATASCAN_1	31-32	DP 7940 _ 9300 mm		0 kPa	
6	1 DC	226	DATASCAN_1	32-33	DP 9300 _ 9800 mm		0 kPa	
7	1 DC	227	DATASCAN_1	33-34	DP 9800 _ 10400 mm		0 kPa	
8	1 DC	228	DATASCAN_1	34-35	DP 10400 _ 11820 mm		0 kPa	
9	1 DC	229	DATASCAN_1	35-36	DP 11820 _ 12410 mm		0 kPa	
0	1 DC	230	DATASCAN_1	36-37	DP 12410 _ 13070 mm		0 kPa	
1	1 DC	231	DATASCAN_1	9300mm	P CYCI ONF INI FT		0 kPa	

Example data sheet:

TIME	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z	32.33
	10	20	30	40	50	60	70	80	90	100	110	120	130	140	150	160	170	180	190	200	210	220	230	240	250	260	270
8.00.00	2.8567	1.8550	0.0016	0.6721	0.2194	0.2027	855.7033	864.0620	2.8819	3.9340	0.0021	3.9361															
8.01.00	2.8653	1.8613	0.0019	0.6727	0.2214	0.1992	855.8229	865.4597	2.9234	3.9459	0.0025	3.9485															
8.02.00	2.8585	1.8544	0.0021	0.6726	0.2208	0.2006	855.8100	866.0333	2.9819	3.9371	0.0028	3.9399															
8.03.00	2.8523	1.8490	0.0000	0.6713	0.2205	0.2017	855.8181	864.9004	2.8704	3.9247	0.0000	3.9247															
8.04.00	2.8678	1.8642	0.0031	0.6719	0.2204	0.2010	855.7743	865.0612	2.8987	3.9455	0.0042	3.9507															
8.05.00	2.8678	1.8655	0.0006	0.6715	0.2217	0.1998	855.8018	865.9154	2.8970	3.9495	0.0008	3.9504															
8.06.00	2.8655	1.8650	0.0024	0.6699	0.2202	0.2003	855.7532	867.0096	2.8683	3.9502	0.0032	3.9534															
8.07.00	2.8435	1.8444	0.0000	0.6702	0.2213	0.1997	855.8019	865.9264	2.8632	3.9161	0.0000	3.9161															
8.08.00	2.8523	1.8495	0.0012	0.6716	0.2196	0.2007	855.8818	864.7214	2.8731	3.9240	0.0016	3.9256															
8.09.00	2.8558	1.8542	0.0006	0.6698	0.2194	0.2014	855.9515	865.8651	2.8797	3.9328	0.0007	3.9360															
8.10.00	2.8670	1.8684	0.0004	0.6732	0.2201	0.2012	856.0325	867.3800	2.9031	3.9587	0.0000	3.9587															
8.11.00	2.8459	1.8408	0.0004	0.6732	0.2204	0.2007	856.0455	867.2815	2.8589	3.9241	0.0006	3.9247															
8.12.00	2.8513	1.8515	0.0000	0.6732	0.2203	0.2002	856.0625	867.2043	2.8763	3.9273	0.0000	3.9273															
8.13.00	2.8543	1.8498	0.0008	0.6734	0.2211	0.1999	856.1686	866.0531	2.8875	3.9315	0.0012	3.9326															
8.14.00	2.8658	1.8605	0.0000	0.6754	0.2209	0.2014	856.2126	865.6756	2.8923	3.9459	0.0000	3.9459															
8.15.00	2.8344	1.8255	0.0032	0.6745	0.2201	0.2032	856.2142	864.1885	2.8397	3.8976	0.0043	3.9018															
8.16.00	2.8904	1.8842	0.0010	0.6746	0.2200	0.1986	856.1840	864.8873	2.9734	3.9770	0.0013	3.9789															
8.17.00	2.8736	1.8689	0.0000	0.6726	0.2196	0.2023	856.1979	866.6602	2.9023	3.9601	0.0000	3.9601															
8.18.00	2.8592	1.8558	0.0055	0.6717	0.2207	0.1999	856.1735	867.4019	2.8833	3.9428	0.0075	3.9503															
8.19.00	2.8680	1.8643	0.0026	0.6731	0.2206	0.2032	856.0876	868.2343	2.8983	3.9578	0.0036	3.9614															
8.20.00	2.8613	1.8596	0.0000	0.6718	0.2201	0.2000	856.0536	867.7916	2.8874	3.9471	0.0000	3.9471															
8.21.00	2.8419	1.8396	0.0015	0.6720	0.2218	0.2025	856.0051	867.1839	2.8586	3.9183	0.0021	3.9204															
8.22.00	2.8630	1.8609	0.0000	0.6721	0.2201	0.2008	856.0325	867.1884	2.8598	3.9474	0.0000	3.9474															
8.23.00	2.8508	1.8473	0.0063	0.6732	0.2203	0.2004	856.1201	867.2991	2.8782	3.9309	0.0086	3.9395															
8.24.00	2.8474	1.8688	0.0000	0.6707	0.2197	0.1995	856.2061	867.3235	2.9023	3.9594	0.0000	3.9594															
8.25.00	2.8614	1.8601	0.0016	0.6705	0.2212	0.1919	856.2321	869.1636	2.8913	3.9155	0.0022	3.9166															
8.26.00	2.8562	1.8562	0.0000	0.6719	0.2203	0.1989	856.1940	868.2880	2.8927	3.9331	0.0000	3.9331															
8.27.00	2.8473	1.8430	0.0114	0.6715	0.2200	0.2015	856.2034	866.2780	2.8642	3.9226	0.0019	3.9245															
8.28.00	2.8465	1.8403	0.0026	0.6735	0.2193	0.1992	856.2564	866.2208	2.8584	3.9212	0.0003	3.9215															
8.29.00	2.8999	1.8637	0.0046	0.6726	0.2195	0.2007	856.2936	866.6605	2.9134	3.9550	0.0062	3.9612															

GC data

Definition: Gas measurement (H2, CO, CH4...)

Format: Excel, 1 sheets with all the data

20190201\GC1_190201_Gasifier_2_1.DATA											
A	B	C	D	E	F	G	H	I	J	K	L
	Acquisition Date & Time	Quantity/H ₂	Quantity/H ₂	Quantity/O ₂	Quantity/N ₂	Quantity/CH ₄	Quantity/CO	Quantity/CO ₂	Quantity/C ₂ H ₄	Quantity/C ₂ H ₆	Quantity/C ₂ H ₂
1	20190201\GC1_190201_Gasifier_2_1.DATA	2019-02-01 11:39:30	1.38	19.59	0.2	16.9	16.54	10.96	15.58	14.2	0.26
2	20190201\GC1_190201_Gasifier_2_2.DATA	2019-02-01 11:43:19	1.02	19.53	0.19	10.71	22.08	9.25	13.29	20.83	0.26
3	20190201\GC1_190201_Gasifier_2_3.DATA	2019-02-01 11:47:09	0	0.02	10.2	72.44	0.11	0.09	0.05	0.01	
4	20190201\GC1_190201_Gasifier_2_4.DATA	2019-02-01 11:50:58	0	0	10.24	72.76	1.06	0.23	0.04		
5	20190201\GC1_190201_Gasifier_2_5.DATA	2019-02-01 11:54:47	0.4	8.12	5.36	40.63	10.91	4.33	6.54	11.01	0.13
6	20190201\GC1_190201_Gasifier_2_6.DATA	2019-02-01 11:58:37	1.11	19.05	0.2	11.25	21.02	9.3	12.52	21	0.27
7	20190201\GC1_190201_Gasifier_2_7.DATA	2019-02-01 12:02:27	1.13	19.86	0.2	11.9	20.9	9.26	12.63	20.08	0.26
8	20190201\GC1_190201_Gasifier_2_8.DATA	2019-02-01 12:06:16	1.18	19.37	0.2	12.24	20.3	10.37	13.22	18.88	0.25
9	20190201\GC1_190201_Gasifier_2_9.DATA	2019-02-01 12:10:06	1.02	19.16	0.19	10.76	22.02	9.88	12.23	21.01	0.25
10	20190201\GC1_190201_Gasifier_2_10.DATA	2019-02-01 12:13:55	1.02	19.62	0.19	10.77	22.09	9.4	12.24	21.17	0.24
11	20190201\GC1_190201_Gasifier_3_1.DATA	2019-02-01 12:17:44	1	20.04	0.19	10.62	20.96	11.56	13.07	19.09	0.23
12	20190201\GC1_190201_Gasifier_3_2.DATA	2019-02-01 12:21:33	1.03	19.8	0.2	10.91	21.37	10.65	12.75	19.87	0.22
13	20190201\GC1_190201_Gasifier_3_3.DATA	2019-02-01 12:25:24	1.07	19.68	0.2	11.74	22.09	8.47	12.39	21.08	0.22
14	20190201\GC1_190201_Gasifier_3_4.DATA	2019-02-01 12:29:13	1.09	19.25	0.19	11.77	23.01	7.14	11.24	22.89	0.27
15	20190201\GC1_190201_Gasifier_3_5.DATA	2019-02-01 12:33:02	0.39	6.44	7.16	53.12	6.79	2.29	3.23	5.39	0.07
16	20190201\GC1_190201_Gasifier_3_6.DATA	2019-02-01 12:36:51	0.01	0.1	10.17	72.38	1.23	0.06	0.05	0.01	
17	20190201\GC1_190201_Gasifier_3_7.DATA	2019-02-01 12:40:40	0.01	0.04	10.24	72.65	0.08	0.04	0.05	0.01	
18	20190201\GC1_190201_Gasifier_3_8.DATA	2019-02-01 12:44:29	0.01	0.02	10.3	73.07	0.07	0.02	0.05	0.02	
19	20190201\GC1_190201_Gasifier_3_9.DATA	2019-02-01 12:48:18	0.01	0.02	10.3	73.14	0.08	0.04	0.05	0.02	
20	20190201\GC1_190201_Gasifier_3_10.DATA	2019-02-01 12:52:09	0.01	0.01	10.3	73.28	0.05	0.04	0.05	0.01	
21	20190201\GC1_190201_Gasifier_4_1.DATA	2019-02-01 12:55:58	0.01	0.01	10.36	73.75	0.04	0.02			
22	20190201\GC1_190201_Gasifier_4_2.DATA	2019-02-01 12:59:47	0.01	0.01	10.35	73.54	0.07	0.04			
23	20190201\GC1_190201_Gasifier_4_3.DATA	2019-02-01 13:03:36	0.01	0.01	10.34	73.47	0.03	0.02			
24	20190201\GC1_190201_Gasifier_4_4.DATA	2019-02-01 13:07:25	0.01	0.01	10.34	73.5	0.02	0.02			
25	20190201\GC1_190201_Gasifier_4_5.DATA	2019-02-01 13:11:15	0.01	0.01	10.36	73.52	0.04	0.03			
26	20190201\GC1_190201_Gasifier_4_6.DATA	2019-02-01 13:15:04	0.01	0.01	10.37	73.57	0.04	0.01			
27	20190201\GC1_190201_Gasifier_4_7.DATA	2019-02-01 13:18:54	0.01	0.01	10.37	73.52	0.01	0.01			
28	20190201\GC1_190201_Gasifier_4_8.DATA	2019-02-01 13:22:43	0.01	0.01	10.37	73.54	0	0.01			
29	20190201\GC1_190201_Gasifier_4_9.DATA	2019-02-01 13:26:32	0.01	0.01	10.35	73.58	0.01	0.02			
30	20190201\GC1_190201_Gasifier_4_10.DATA	2019-02-01 13:30:21	0.01	0.01	10.36	73.55	0.02	0.02			
31	20190201\GC1_190201_Gasifier_5_1.DATA	2019-02-01 13:34:10	0.01	0.01	10.34	73.54	0.05	0.02			

SPA data

Definition: Aromatics and other heavy molecules

Format: Excel, multiple sheets with all the data

	A	B	C	D	E	F	G	H	I	J
1	Index	Name	Time	Quantity	Height	Area	Area %			
2			[Min]	[mg/l]	[uV]	[uV.Min]	[%]			
3	1	Benzene	1.5	6.48	933.6	42.6	5.568			
4	2	UNKNOWN	1.98	3.61	226.5	24	3.142			
5	3	Toluene	2.29	2.4	273.1	16	2.087			
6	4	Styrene	4.87	1.49	184.7	10.1	1.318			
7	5	Indene	8.55	1.12	205.5	7.7	1.006			
8	6	Naphthalene	11.68	8.01	1884.9	56.9	7.436			
9	7	2-MethylNaphthalene	13.72	0.87	167	6.2	0.809			
10	8	Intern standard 1	14.02	90.68	12621.1	391.8	51.195			
11	9	1-MethylNaphthalene	14.28	0.96	178.6	6.8	0.891			
12	10	Biphenyl	15.9	0.92	160.1	6.7	0.87			
13	11	UNKNOWN	17.12	0.92	135.1	6.1	0.798			
14	12	Acenaphthylene	18.29	3.71	561.1	24.5	3.201			
15	13	Fluorene	21.2	0.9	150.8	6.6	0.859			
16	14	Phenantrene	25.14	2.61	612.8	19.2	2.511			
17	15	Anthracene	25.25	0.62	139.4	4.7	0.616			
18	16	UNKNOWN	29.13	0.98	205.6	6.6	0.856			
19	17	Fluoranthene	29.21	0.09	19.5	0.7	0.093			
20	18	Pyrene	30	1.49	339.6	11.4	1.485			
21	19	UNKNOWN	31.01	1	162.4	6.6	0.867			
22	20	UNKNOWN	32.61	0.89	117	5.9	0.776			
23	21	UNKNOWN	34.42	15.64	3385.2	104.2	13.615			
24										
25	Total	-	-	145.41	22663.8	765.3	100			
26										
27										
28										
29										
30										
31										
32										

430_190206_190201_G14_F_44_2430_190206_190201_G14_F_44_1430_190206_190201_G14_C_43_3

Any file needed for research (Excel, photos, etc.). The program will only store but not analyse.

Point Window

The *Point window* includes the main information of the Point: name, start and end date and time and comments. This information can be changed if needed in *Modify info*.

It also includes a table Data already linked to the point (*Data linked*) Including multiple actions.

If more data wants to be included/linked to the point, select the type of data and click *Link Data*.

Season 1/ Experiment 1/Point 1

NamePoint0

Commentsthis is the point 0

MODIFY INFO

Status: Ready!

Data Linked

Index	Type	StartDate	EndDate	Delay	Ientrie
-------	------	-----------	---------	-------	---------

VIEW DATA

ANALYZE DATA

DELETE DATA

Date (YYYY-MM-DD)Time (HH:MM:SS)

Date start2020-10-18:

Date end2020-10-117:00

Link data to point

Select type

AUTOMATIC
SCADA
GC1
INFERNO
GPA

LINK DATA

OK

Action 1: Modify info

In order to modify the information of a Point (names and comments) click *Modify Info*. The text will turn white (same as the [Action 1: Modify Experiment Attributes](#)) and you can modify any attribute, when finish, click *Modify Info* again.

Action 2: View Data

The table Data linked e shows the existing data linked to the point. If you want to see the raw data click *View Data*.

Action 3: Delete Data

The table Data linked shows the existing data linked to the point. If you want to delete the raw data click *Delete Data*.

Action 4: Analyze data

The table Data linked shows the existing data linked to the point. If you want to do calculation or analysis on this data click *Analyse Data*.

This function is under construction



Action 5: Link data

If more data wants to be included/linked to the point, select the type of data and click *Link Data*.

For all data types is the same procedure, write down the time when the points take place (normally a point has a duration from 30 min to 1 hour and has stable operation conditions).

If there is a delay of the GC computer (not connected to internet) with respect the SCADA computer must be added (For example if GC time is 10:05 and SCADA 10:02, the delay is +3).

MainWindow

Read Data

Collect Data

	Date (YYYY-MM-DD)	Time (HH:MM:SS)
Date start	<input type="text" value="2019-02-01"/>	<input type="text" value="11:55:00"/>
Date end	<input type="text" value="2019-02-01"/>	<input type="text" value="12:27:00"/>
Time Delay resp. SCADA (HH:MM:SS)	<input type="text" value="00:03:00"/>	

Delay of data selected with regards SCADA. For example,
if GC computer time is 10:05 and KC computer time is 10:02,
the delay is of +3min

3. How to implement code changes

Please refer to: \Project_KC\docs\build\html and open file *index* in chrome or Firefox.